


# IntelliScan Pulse Oximeter Model CVD-M1002 & MBT1003

## IntelliScan Pulse Oximeter, Model CVD-M1002:

- IntelliScan Oxygen Saturation (SpO<sub>2</sub>) Monitor
- Pulse rate monitor
- Bluetooth technology for data transmission from oximeter to mobile phone (email/Text/Skype)
- Dual direction display
- Auto power shut-off
- Battery level display indicator
- SpO<sub>2</sub>  
Measurement range : 35 ~ 100%  
Accuracy : ± 2% (80 ~ 100%)
- Pulse rate:  
Measurement range : 25 ~ 250bpm  
Accuracy : ± 2bpm
- Display : OLED LCD
- Power : 2 x AAA 1.5V alkaline batteries
- Working time: Continuously 30 hours Approx.
- Accessories : Soft Bag & User Manual packaged in Box.
- Dimension : 58x34x30mm (2.3"x1.3"x1.2") Approx.
- Weight : 50g (2oz) Approx.



IntelliScan Pulse Oximeter  
Model CVD-M1002



### Disclaimer:

For general wellness and health applications including sports, fitness, and relaxation management. Not intended to be used in prevention, diagnosis, or treatment of any disease. The product is designed to provide instant information to consumer or organization for appropriate action to reach out a licensed medical practitioner. Competent medical advice should be sought if there are concerns around illness, and further checks may be required under the supervision of a competent medical practitioner. You should not use this information as self-diagnosis or for treating a health problem or disease.

# IntelliScan Pulse Oximeter Model CVD-M1002 & MBT1003

## IntelliScan BT Pulse Oximeter, CVD-BT1003:

- IntelliScan Oxygen Saturation (SpO<sub>2</sub>) Monitor
- Pulse rate monitor
- Bluetooth technology for data transmission from oximeter to mobile phone (email/Text/Skype)
- Dual direction display
- Auto power shut-off
- Battery level display indicator
- SpO<sub>2</sub>:  
Measurement range : 35 ~ 100%  
Accuracy : ± 2% (80 ~ 100%)
- Pulse rate:  
Measurement range : 25 ~ 250bpm  
Accuracy : ± 2bpm
- Display : OLED LCD
- Power : 2 x AAA 1.5V alkaline batteries
- Working time : 15 hours continuously
- Accessories : Soft Bag & User Manual packaged in Box.
- APP : IntelliScan OxyCare (Android/iOS)
- Dimension : 58x34x30mm (2.3"x1.3"x1.2") Approx.
- Weight : 50g (2oz) Approx.



### General Notes:

1. IntelliScan Pulse Oximeter, Model CVD-M1002 is designed and manufactured IAW general engineering practices, USFDA CFR's, and user friendly features.
2. The medical and healthcare products are the backbone of our society and therefore our devices offer the highest level of performance and comply with rigorous requirements.
3. Intelliscan Systems, associated and support products are designed and manufactured in accordance with the intent of specifications such as MIL-STD-810, MIL-I-45208, MIL-STD-883C, MIL-STD-45662, MIL-STD-129, MIL-STD-130, MIL-STD-105, MIL-STD-454, MIL-STD-595 and other applicable specifications such as SAE, ASTM and as required by a contract.

### Disclaimer:

For general wellness and health applications including sports, fitness, and relaxation management. Not intended to be used in prevention, diagnosis, or treatment of any disease. The product is designed to provide instant information to consumer or organization for appropriate action to reach out a licensed medical practitioner. Competent medical advice should be sought if there are concerns around illness, and further checks may be required under the supervision of a competent medical practitioner. You should not use this information as self-diagnosis or for treating a health problem or disease.